



Data Sheet

Product Name: Human Liver Tissue (Non- cancer) FFPE Sections
SKU#: HuPS-08007A

Product type: FFPE Sections

Species: Human

Tissue Name: Liver

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

Tissue	Pathology Diagnosis	Gender/Age	IHC Data
Human Liver	Benign liver tissue with mild portal fibrosis, mild portal and lobular chronic inflammation.	Female/43	NA

Package size

- 1 section per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.