



## Data Sheet

**Product Name:** Human Liver Metastatic Adenocarcinoma FFPE Sections

**SKU#:** HuPS-08015T

**Product type:** FFPE Sections

**Species:** Human

**Tissue Name:** Liver

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5  $\mu$ m in thickness and mounted each on positively charged glass slides.

**Quality Control:**

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

**Donor and Histopathology Information**

Tissue	Pathology Diagnosis	Gender/Age (year)	% Tumor area	Tumor Grade	TMN Stage	IHC Data
Human liver	Metastatic moderately differentiated adenocarcinoma, Benign liver tissue with mild portal fibrosis, mild portal chronic inflammation	Male/68	55%	II	N/A	

**Package size**

- 1 section per slide
- 5 slides per package

**Applications**

Histology stain, IHC/IF, FISH and ISH

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.