



**Data Sheet**

**Product Name:** Human Sinus (normal) FFPE Sections

**SKU#:** HuPS-15001A

**Product type:** FFPE Sections

**Species:** Human

**Tissue Name:** Normal Sinus

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

**Quality Control:**

- 1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
- 2. IHC tested: yes

**Donor and Histopathology Information**

Tissue	Pathology Diagnosis	Gender/Age (year)	IHC Data
Human Sinus	Normal Sinus	Male/67	

**Package size**

- 1 section per slide
- 5 slides per package

**Applications**

Histology stain, IHC/IF, FISH and ISH

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.