



### Data Sheet

Product Name: Human Breast Cancer (ER +/-PR-/Her2 -), FFPE Sections  
SKU#: HuPS-02037T1

**Product type:** FFPE Sections

**Species:** Human

**Tissue Name:** Breast

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

**Quality Control:**

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

**Donor and Histopathology Information**

Tissue	Pathology Diagnosis	Gender/Age (year)	Tissue of Origin	% Tumor area	Tumor Grade	TMN Stage	IHC Data
Breast Cancer	Human breast Invasive ductal carcinoma	Female/47	Human breast, right	40%	III	N/A	ER (+/-, 5%) PR(-) Her2 (+/- Borderline) Ki67(70%)

**Package size**

- 1 section per slide
- 5 slides per package

**Applications**

Histology stain, IHC/IF, FISH and ISH

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.