



Data Sheet

Product Name: PD-L1 Positive (Human Lung Adenocarcinoma) FFPE Sections

SKU#: HuPS-03001T-1

Product type: FFPE Sections

Species: Human

Tissue Name: Lung

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

Organ& Tissue	Pathology Diagnosis	Gender/Age	Smoking	TMN/Stage	Grade	Biomarker
Lung Cancer	Adenocarcinoma	Female/73			2	PD-L1 (++), Ki67(+) 60%, CK7 (+), EGFR(+), TTF-1 (+)

Package size

- 1 section per slide
- 5 slides per package

Applications

PD-L1 IHC/IF positive control

Storage

2-8°C

Note: This product is intended for research purposes only.