



## Data Sheet

**Product Name:** Human Pancreas Tissue (Normal, Ampulla of Vater) FFPE Sections

**SKU#:** HuPS-12003A

**Product type:** FFPE Sections

**Species:** Human

**Tissue Name:** Pancreas

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

**Quality Control:**

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

**Donor and Histopathology Information**

| Tissue   | Pathology Diagnosis                                     | Gender/Age (year) | Tissue of Origin                   | IHC Data |
|----------|---|-------------------|------------------------------------|----------|
| Pancreas | Normal architecture and cellularity of ampulla of vater | Male/74           | Human Pancreas of ampulla of vater |          |

**Package size**

- 1 section per slide
- 5 slides per package

**Applications**

Histology stain, IHC/IF, FISH and ISH

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.