



## Data Sheet

**Product Name:** Positive Control Slides of Human Alzheimer's Disease Brain FFPE Sections for Tau IHC/IF and Thioflavin S Staining

**SKU#:** HuPS-AD+

**Product type:** FFPE Sections

**Species:** Human

**Tissue Name:** Brain

**Slide Preparation:** Human brain tissue from an Alzheimer's Disease autopsy was immediately fixed in formalin following excision and subsequently embedded in paraffin. The tissue sections, each 5  $\mu$ m thick, were mounted onto positively charged glass slides.

**Quality Control:** At least one of the tissue slides from each lot was stained with Tau IHC and Thioflavin S staining to ensure the quality

**Package size**

1 section per slide  
5 slides per package

**Applications**

Tau IHC and Thioflavin S staining positive control

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.