



### Data Sheet

**Product Name:** Human Ovary High Grade Serous Carcinoma (HGSC) FFPE Sections

**SKU#:** HuPS-14016T1

**Product type:** FFPE Sections

**Species:** Human

**Tissue Name:** Ovary

**Slide Preparation:** Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

**Quality Control:**

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

**Donor Information:**

Organ& Tissue	Pathology Diagnosis	Gender/Age	% Tumor Area	Grade	TMN Stage	Biomarkers (IHC)
Human Ovary	Human ovary high grade serous carcinoma (HGSC)	Female/67	80	III	NA	ER(+),Pax-8 (+),WT-1(+),p16(+),p53(+),Ck7 (+), CK20(-)

**Package size**

- 1 section per slide
- 5 slides per package

**Applications**

Histology stain, IHC/IF, FISH, ISH, RNAscope ISH and DNA/RNA isolation for NGS

**Storage**

2-8°C

**Note:** This product is intended for research purposes only.