



Data Sheet

Product Name: Human Ovarian Clear-Cell Carcinoma FFPE Sections

SKU#: HuPS-14018T

Product type: FFPE Sections

Species: Human

Tissue Name: Ovary

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

Organ& Tissue	Pathology Diagnosis	Gender/Age	% Tumor Area	Grade	TMN Stage	Biomarkers (IHC)
Human Ovary	Human ovarian clear-cell carcinoma or ovarian clear-cell adenocarcinoma	Female/67	80	II	pT1aNxMx	ER(-),PR(-),Pax8(+),WT1(-),p53(-),Ki67(70%)

Package size

- 1 section per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH, ISH, RNAscope ISH and DNA/RNA isolation for NGS

Storage

2-8°C

Note: This product is intended for research purposes only.