



Data Sheet

Product Name: Human Cervical Adenocarcinoma FFPE Sections

SKU#: HuPS-16003T1

Product type: FFPE Sections

Species: Human

Tissue Name: Cervix

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor Information:

Organ& Tissue	Pathology Diagnosis	Gender/Age	% Tumor Area	Grade	TMN Stage	Biomarkers
Human Cervix	Cervical Adenocarcinoma with mucinous differentiation	Female/38	82%	II	NA	P16(+),P40(-),CK5(+), PD-L1 (+)

Package size

- 1 section per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH, ISH, RNAscope ISH and DNA/RNA isolation for NGS

Storage

2-8°C

Note: This product is intended for research purposes only.