



Data Sheet

Product Name: Human Airway & Blood Vessel FFEP Sections for Positive Control of PAS, Collagen or Elastic Fiber Staining

SKU#: HuPS-Airway

Product type: FFPE Sections

Species: Human

Slide Preparation: Human lung airway and blood vessel tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 μ m in thickness and mounted each on positively charged glass slides.

Quality Control:

At least one of the tissue slides from each lot was stained with PAS to ensure the quality.

Package size

1 section per slide
5 slides per package

Applications

PAS, Manson's trichrome staining and collagen or elastic fiber staining positive control.

Storage

2-8°C

Note: This product is intended for research purposes only.