



Data Sheet

Product Name: Positive Control FFPE Slides for *In Situ* Fungus Detection

SKU#: MsPs-718Fungus+

Product type: FFPE Sections

Species: Mouse

Slide Preparation: Lung tissue from FVB mice infected with fungus was fixed in formalin for 24 hours immediately after excision. The tissue was then processed, embedded in paraffin, and sectioned at a thickness of 5 μ m. Each section was mounted on a positively charged glass slide.

Quality Control:

At least one of the tissue slides from each lot was stained with Grocott-Gomori's Methenamine Silver (GMS) Staining and PAS to ensure the quality.

Package size

1 section per slide
5 slides per package

Applications

Positive Control FFPE Slides for *In Situ* Fungus Detection

Storage

2-8°C

Note: This product is intended for research purposes only.