



Data Sheet

Product Name: Human Ovary B-Lymphoblastic leukemia (ALL) FFPE Sections

SKU#: HuPS-18007T

Product type: FFPE Sections

Species: Human

Tissue Name: Ovary

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 μ m in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor and Histopathology Information

Tissue	Pathology Diagnosis	Gender/Age (year)	% Tumor area	Tumor Grade	TMN Stage	IHC Data
Ovary	Human Ovary B-Lymphoblastic leukemia (ALL)	Female/19	90%	N/A	N/A	Diffusely (+): PAX-5, CD10, TdT Weak (+): CD79, CD34 and CD45 Neg: CD20, CD3 and MPO

Package size

- 1 section per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.