



Data Sheet

Product Name: Human Prostate Adenocarcinoma FFPE Sections
SKU#: HuPS-25001T

Product type: FFPE Sections

Species: Human

Tissue Name: Prostate

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor and Histopathology Information

Tissue	Pathology Diagnosis	Gender/Age (year)	Tissue of Origin	% Tumor area	Tumor Grade	TMN Stage	IHC Data
Prostate Cancer	Human Prostate Adenocarcinoma, Acinar type, Gleason score 9 (4+5)	Male/73	Human right posterior prostate	50%	II	pT3b, PN1 Mx	N/A

Package size

- 1 section per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.