



Data Sheet

Product Name: Human Prostate Adenocarcinoma FFPE Sections
SKU#: HuPS-25002T

Product type: FFPE Sections

Species: Human

Tissue Name: Prostate

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: NA

Donor and Histopathology Information

| Tissue | Pathology Diagnosis | Gender/Age (year) | Tissue of Origin | % Tumor area | Tumor Grade | TMN Stage | IHC Data |
|-----------------|---|-------------------|--------------------------------|--------------|-------------|--------------|----------|
| Prostate Cancer | Human Prostate Adenocarcinoma, Acinar type, Gleason score 7 (4+3) | Male/66 | Human right posterior prostate | 31 | II | pT3a, PNx Mx | NA |

Package size

- 1 section per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.