



Data Sheet

Product Name: Human Diffuse Large B-Cell Lymphoma (DLBCL) FFPE Sections

SKU#: HuPS-18003T

Product type: FFPE Sections

Species: Human

Tissue Name: Lymph Node

Slide Preparation: Human tissue was fixed in formalin immediately after excision and embedded in paraffin. The tissue sections were 5 µm in thickness and mounted each on positively charged glass slides.

Quality Control:

1. At least one of the tissue slides from each lot was stained with H & E to ensure the quality
2. IHC tested: yes

Donor and Histopathology Information

Tissue	Pathology Diagnosis	Sample Procedure	Treatment Status	Gender/Age (year)	Tumor Grade	TMN Stage	IHC Data
Human lymph node	Diffuse Large B-Cell Lymphoma (DLBCL)	Surgical Resection	Treatment Naïve	Male/52	N/A	N/A	CD20+, CD3-, MIB Index 65-70%

Package size

- 1 section per slide
- 5 slides per package

Applications

Histology stain, IHC/IF, FISH and ISH

Storage

2-8°C

Note: This product is intended for research purposes only.